

Mölnlycke[®]
Exufiber[®] Ag+







STERILE EO

CE 2797

Manufacturer

 Mölnlycke Health Care AB
Gamlestadvägen 3C, Box 13080
SE-402 52 Göteborg, Sweden

www.molnlycke.com

en For low exuding wounds

For moderately exuding wounds

For highly exuding wounds

no Brukes på lett væskende sår

Brukes på moderat væskende sår

Brukes på sterkt væskende sår

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Exufiber® Ag+

Gelling Fibre Dressing with silver

en

Product description

Exufiber Ag+ gelling fibre dressing with silver is a sterile, soft, nonwoven dressing made from highly absorbent and gel forming polyvinyl alcohol (PVA) fibres, coated with silver sulphate on both sides. Exufiber Ag+ has high absorption and retention capacity of wound exudate and maintains a moist wound environment. The structural integrity of Exufiber Ag+ supports removal in one piece.

In contact with fluid Exufiber Ag+ releases silver ions to the wound bed and provides an antimicrobial action. It inactivates a wide range of wound related pathogens, such as bacteria, fungi and mould, as shown *in vitro*. By reducing the number of microorganisms, Exufiber Ag+ can prevent the formation of and reduce bacterial biofilm, as shown *in vitro* and *in vivo*. Exufiber Ag+ may be used as an effective barrier to microbial penetration of the dressing.

Intended use

Exufiber Ag+ is intended to be used in the following medium to high exuding wounds:

- Venous leg ulcers
- Diabetic foot ulcers

Indicated wear time: up to seven (7) days.

Instructions for use

Shallow wounds

Application:

Cleanse the wound with saline solution or water according to clinical practice.

1. Dry the surrounding skin thoroughly.
2. Choose the correct size of dressing to be able to cover the entire wound. For best result, Exufiber Ag+ should overlap the dry surrounding skin by at least 1 – 2cm for the smaller sizes (up to 10x10cm) and 5cm for larger sizes. The dressing will swell, causing contraction of the dressing edges, as it absorbs wound fluid and starts gelling.
3. Apply a dry Exufiber Ag+ dressing to the wound.
4. Fixate with an appropriate secondary dressing. The choice depends on the exudate level.
5. Compression therapy may be used in conjunction with Exufiber Ag+.

Deep wounds and cavities

Application:

Cleanse the wound with saline solution or water according to clinical practice.

1. Dry the surrounding skin thoroughly.
2. Loosely pack ribbon or pad into the wound to allow room for swelling of the dressing. When using the ribbon dressing, cut to appropriate length leaving a small overhang of 2– 3cm outside the wound for easy retrieval.
3. Cover with an appropriate secondary dressing. The choice depends on the exudate level.
4. Fixate (or cover) with an adhesive dressing or a bandage as appropriate.

Dressing change and removal

The wound and Exufiber Ag+ should be inspected regularly. Exufiber Ag+ can be left in place for up to 7 days, depending on wound condition or as indicated by clinical practice.

The secondary dressing should be changed based on exudate level or before becoming fully saturated.

1. Remove Exufiber Ag+ by gently cleansing/flushing with saline solution or water according to clinical practice. Any non-gelled material will moisten in contact with the saline.
2. If Exufiber Ag+ dries out and is difficult to remove, it should be moistened with saline solution or water according to local policies and allowed to soak until it lifts easily. It may take several minutes for Exufiber Ag+ to transform into a gel. Remove the dressing by gently cleansing/flushing.

Exufiber Ag+ is intended for short-term use up to 4 weeks. For long-term use, a clinical assessment by a physician is recommended

Contraindications

- For external use only.

Precautions

- In case of signs of clinical infection, consult a health care professional. Exufiber Ag+ does not replace the need for systemic or other adequate infection treatment.
- Exufiber Ag+ should be used under the supervision of a qualified health care professional.
- Do not use on patients with a known sensitivity to silver or any other contents of the dressing.
- Clinicians/ Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Exufiber Ag+ may cause transient discoloration of the wound bed and surrounding skin.
- Prior to commencing radiation therapy, remove Exufiber Ag+ if product is present in the treatment area. A new dressing can be applied following treatment.

en

- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).
- Do not use Exufiber Ag+ together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- Other than saline solution or water, the interaction of cleansing agents in combination with Exufiber Ag+ has not been demonstrated.
- The interaction of Exufiber Ag+ with topical treatments has not been demonstrated.

Storage and disposal

Exufiber Ag+ should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any colour variances in Exufiber Ag+ do not affect product performance or safety.

Disposal should be handled according to local environmental procedures.

Other information

Exufiber Ag+ contains polyvinyl alcohol fibres, hydroxypropylcellulose and 0.2mg/cm² silver.

The PVA fibres used in the product may change colour when it is exposed to light, air and/or heat. The colour change has no influence on product properties when used before expiry date.

Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.

If the product is used after the expiry date product properties cannot be ensured.

Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.

Exufiber is a registered trademark of Mölnlycke Health Care AB.

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CANADA Mölnlycke Health Care Inc.,
2010 Winston Park Drive, Suite 100,
Oakville, Ontario, L6H 5R7
Tel: +1 905 829 1502

CHINA Mölnlycke Health Care
(Shanghai) Co., Ltd.
RM 629, No.8 Hua Jing Road,
Wai Gao Qiao FTZ Shanghai, China
Tel: +86 10 5128 8571

CZECH REPUBLIC Mölnlycke Health Care s.r.o.
Hajkova 2747/22, 130 00 Prague 3
Tel: +420 221 890 517 (Reception)
Tel: +420 221 890 511 (Customer Service)

DENMARK Mölnlycke Health Care ApS,
Gydevang 33, DK-3450 Allerød
Tel: +45 48 16 82 68
Customer Service Tel: +80 886 810

ESTONIA Mölnlycke Health Care OÜ,
Vabaõhumuuseumitee 4-67, EE-13522 Tallin
Tel: +372 671 1520

FINLAND Mölnlycke Health Care Oy,
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki
Tel: +358 201 622 300
Customer Service Tel: 0800 113 442

FRANCE Mölnlycke Health Care SAS,
13 Allée du Château Blanc
59290 WASQUEHAL
Tel: +33 320 122 555
Customer Service Tel: 0800 910 292

GERMANY Mölnlycke Health Care GmbH
Grafenberger Allee 297
DE-40237 Düsseldorf
Tel: +49 211 92 08 80
Customer Service Tel: 0800 186 21 80

HUNGARY Mölnlycke Health Care Kft
HU-1134 Budapest, Dévai u. 26-28. III. em.
Tel: +36 1 477 3080

INDIA Mölnlycke Health Care India Pvt Ltd
Universal Business Centre
Second Floor, Universal Trade Tower
Sector 49, Gurgaon Sohna Road
Gurgaon 122 018, Haryana
Tel: +91 124 4696 008

ITALY Mölnlycke Health Care s.r.l.,
Via Marsala 40/C, IT-21013 Gallarate (VA)
Tel: +39 0331 714 011

JAPAN Mölnlycke Health Care KK,
4F Conceria Nishi-Shinjuku Tower's West
6-20-7 Nishi-Shinjuku,
Shinjuku-Ku, Tokyo 160-0023
Tel: +81 3 6914 5004

KOREA Mölnlycke Health Care Korea Co., Ltd
2nd floor, Hansol Richville.
46 Wiryeseong-daero,
Songpa-gu, Seoul, 05627
South Korea
Tel: +82 2 3789 1402

LATVIA Mölnlycke Health Care AB,
Dzelzavas 120 M, Riga LV-1021
Tel: +372 6711 522

LITHUANIA Mölnlycke Health Care AB,
J. Galvydzio g. 3, LT 08236 Vilnius
Tel: +370 5 2274 58 37

NETHERLANDS Mölnlycke Health Care B.V.
Postbus 3196, 4800 DD Breda
Tel: +31 76 521 9663

NORWAY Mölnlycke Health Care AS
Postboks 6229 Etterstad, NO-0603 Oslo
Tel: +47 22 70 63 70
Customer Service Tel: 800 161 78

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Mölnlycke Health Care Ltd.
USA Mölnlycke Health Care US, LLC