
Mölnlycke®

Mepilex® Transfer 

WITH SAFETAC® TECHNOLOGY

Manufacturer

 Mölnlycke Health Care AB
Gamlestadvägen 3C, Box 13080
SE-402 52 Göteborg, Sweden

www.molnlycke.com

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Mepilex® Transfer Ag with Safetac® technology

Antimicrobial soft silicone exudate transfer dressing

Product description

Mepilex Transfer Ag is a soft silicone wound contact layer that absorbs and transfers exudate, maintains a moist wound healing environment and has antimicrobial properties.

Mepilex Transfer Ag contains silver sulphate which when in contact with fluid will release silver ions to create an effective bacterial barrier and inactivate a wide range of wound related pathogens (bacteria and fungi), as shown in vitro. By reducing the number of micro-organisms, Mepilex Transfer Ag may also reduce odour.

Mepilex Transfer Ag has been shown to inactivate wound related pathogens up to 14 days in vitro.

Mepilex Transfer Ag consists of:

1. a Safetac® adhesive layer
2. a compressed polyurethane foam containing silver sulphate and activated carbon as colouring agent

Safetac® is a unique adhesive technology that minimises pain to patients and trauma to intact skin or wounds.

Mepilex Transfer Ag is soft and conformable making it easy to keep the dressing in contact with the wound surface and the surrounding skin, even in awkward or uneven areas.

Mepilex Transfer Ag maintains a moist wound environment in combination with an appropriate secondary dressing.

Mepilex Transfer Ag can be cut to suit various wound shapes and locations.

Dressing material content:

Polyurethane, silicone, polyolefin, polyacrylate, silver sulphate

Intended use

Mepilex Transfer Ag is designed for a wide range of low to high exuding and difficult to dress leg and foot ulcers, pressure ulcers, malignant wounds, partial thickness burns and donor sites, where there is a risk of infection.

Mepilex Transfer Ag may be used on infected wounds as part of a treatment regimen under supervision of a qualified health care professional.

Mepilex Transfer Ag can be used under compression bandaging.

Instructions for use

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Remove product from sterile pack using aseptic technique. Cut to appropriate size (if needed). Remove the release films and apply the product with the adherent side to the wound. Do not stretch.
4. For best adherence to surrounding skin cut Mepilex Transfer Ag to a size that overlaps the wound area by at least 2 cm for the smaller product sizes and 5 cm for the larger product sizes.
If more than one piece of Mepilex Transfer Ag is required overlap the edges.
5. Apply a secondary dressing appropriate for the level of exudate. The secondary dressing should exceed the edges of Mepilex Transfer Ag. Fixate to secure close contact with the wound.

Mepilex Transfer Ag is intended for short-term use up to 4 weeks. For continued use, a re-assessment by a physician is recommended.

Frequency of change

Mepilex Transfer Ag may be left in place for up to 14 days depending on the patient, condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

Precautions

- Mepilex Transfer Ag should be used under the supervision of a qualified health care professional.
- Do not use on patients with a known sensitivity to silver or any other contents of the dressing.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Mepilex Transfer Ag may cause transient discoloration of the wound bed and surrounding skin.
- In the event of clinical infection Mepilex Transfer Ag does not replace the need for systemic therapy or other adequate infection treatment.
- Do not use Mepilex Transfer Ag during radiation treatment or examinations e.g. X-ray, ultrasound, diathermy or Magnetic Resonance Imaging.

- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).
- Do not use Mepilex Transfer Ag together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- Other than saline solution or water, the interaction of cleansing agents in combination with Mepilex Transfer Ag has not been demonstrated.
- The interaction of Mepilex Transfer Ag with topical treatments has not been demonstrated.
- For external use only.
- Do not reuse. If reused, performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.
- Do not use after expiry date. If the product is used after the expiry date product properties cannot be ensured.

Storage and disposal

- Mepilex Transfer Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any colour variances in Mepilex Transfer Ag do not affect product performance or safety.
- Disposal should be handled according to local environmental procedures.

Silver content information

Mepilex Transfer Ag contains 1.2 mg/cm² silver.

Mepilex® and Safetac® are registered trademarks of Mölnlycke Health Care AB

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CANADA Mölnlycke Health Care Inc.,
2010 Winston Park Drive, Suite 100,
Oakville, Ontario, L6H 5R7
Tel: +1 905 829 1502

CHINA Mölnlycke Health Care
(Shanghai) Co., Ltd.
RM 629, No.8 Hua Jing Road,
Wai Gao Qiao FTZ Shanghai, China
Tel: +86 10 5128 8571

CZECH REPUBLIC Mölnlycke Health Care s.r.o.
Hajkova 2747/22, 130 00 Prague 3
Tel: +420 221 890 517 (Reception)
Tel: +420 221 890 511 (Customer Service)

DENMARK Mölnlycke Health Care ApS,
Gydevang 33, DK-3450 Allerød
Tel: +45 48 16 82 68
Customer Service Tel: +80 886 810

ESTONIA Mölnlycke Health Care OÜ,
Vabaõhumuuseumitee 4-67, EE-13522 Tallin
Tel: +372 671 1520

FINLAND Mölnlycke Health Care Oy,
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki
Tel: +358 201 622 300
Customer Service Tel: 0800 113 442

FRANCE Mölnlycke Health Care SAS,
13 Allée du Château Blanc
59290 WASQUEHAL
Tel: +33 320 122 555
Customer Service Tel: 0800 910 292

GERMANY Mölnlycke Health Care GmbH
Grafenberger Allee 297
DE-40237 Düsseldorf
Tel: +49 211 92 08 80
Customer Service Tel: 0800 186 21 80

HUNGARY Mölnlycke Health Care Kft
HU-1134 Budapest, Dévai u. 26-28. III. em.
Tel: +36 1 477 3080

INDIA Mölnlycke Health Care India Pvt Ltd
Universal Business Centre
Second Floor, Universal Trade Tower
Sector 49, Gurgaon Sohna Road
Gurgaon 122 018, Haryana
Tel: +91 124 4696 008

ITALY Mölnlycke Health Care s.r.l.,
Via Marsala 40/C, IT-21013 Gallarate (VA)
Tel: +39 0331 714 011

JAPAN Mölnlycke Health Care KK,
4F Conciera Nishi-Shinjuku Tower's West
6-20-7 Nishi-Shinjuku,
Shinjuku-Ku, Tokyo 160-0023
Tel: +81 3 6914 5004

KOREA Mölnlycke Health Care Korea Co., Ltd
2nd floor, Hansol Richville.
46 Wiryeseong-daero,
Songpa-gu, Seoul, 05627
South Korea
Tel: +82 2 3789 1402

LATVIA Mölnlycke Health Care AB,
Dzelzavas 120 M, Riga LV-1021
Tel: +372 6711 522

LITHUANIA Mölnlycke Health Care AB,
J. Galvydzio g. 3, LT 08236 Vilnius
Tel: +370 5 2274 58 37

NETHERLANDS Mölnlycke Health Care B.V.
Postbus 3196, 4800 DD Breda
Tel: +31 76 521 9663

NORWAY Mölnlycke Health Care AS
Postboks 6229 Etterstad, NO-0603 Oslo
Tel: +47 22 70 63 70
Customer Service Tel: 800 161 78

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Mölnlycke Health Care Ltd.
USA Mölnlycke Health Care US, LLC