

CAD Ref: 10530
 Dimensions: 69.5 x 63.5 x 57mm



Mölnlycke®

Ref no: 47466-11

Component Type: Retail box

Created by: Natalie Whiteside

Refer to ARS PD-527605 and Master Text IFU PD-571437

Mölnlycke Health Care, Gamlestadsvägen 3C, SE-402 52, Göteborg.

Colour information

PMS 361

Varnish free area 37mm x 50mm

PMS 355

PMS Cool Grey 10

Cutter



Two way stretch tubular bandage
 For dressing retention, wet and dry wrapping. Select the correct size and cut to length, open the bandage and stretch over the affected area.
Precautions: Keep clean. If soiled with bodily fluids, discard appropriately and replace. Avoid naked flames and ignition sources, especially when used with paraffin based products. If condition worsens discontinue use and seek medical advice. Single-patient, multi use.
 If any serious incident has occurred in relation to the use of Tubifast, it should be reported to Mölnlycke Health Care and to your local competent authority.
Fibre content: Viscose, Nylon, Elastane.
 Mölnlycke Health Care AB
 Gamlestadsvägen 3C, Box 13080, SE-402 52 Göteborg, Sweden
 The Mölnlycke and Tubifast trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 www.molnlycke.com

GREEN LINE 3m

5 055158 009042

NOMINAL LIMB SIZE

Red Line	9-18cm
Green Line	14-24cm
Blue Line	24-40cm
Yellow Line	35-64cm
Purple Line	64-130cm

GREEN LINE 3m

Mölnlycke®
 Tubifast®

Two way stretch tubular bandage

Australian Sponsor Address:
 Mölnlycke Health Care Pty. Ltd.
 12 Narabang Way
 Betrose, NSW 2085
 Australia

REF 992485

GREEN LINE 3m

Mölnlycke®
 Tubifast®

Two way stretch tubular bandage

1 nominal 5cm x 3m roll

