


Mölnlycke® dressing cutting guide



Mepilex®



Tip of toe
10x10 cm





Hand
10x10 cm
or 10x20 cm



Around tubes
10x10 cm
or 10x20 cm



Elbow
10x10 cm



Ear
10x10 cm



Between toes
10x10 cm



Heel
10x10 cm



General all purpose
for younger patients
10x10 cm





See all video tutorials at:
<https://Lead.me/bcKOMZ>
Or scan code above

Proving it every day

At Mölnlycke®, we deliver innovative solutions for managing wounds, improving surgical safety and efficiency and preventing pressure ulcers. Solutions that help achieve better outcomes and are backed by clinical and health-economics evidence.

In everything we do, we are guided by a single purpose: to help healthcare professionals perform at their best. And we're committed to proving it every day.



Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00. The Mölnlycke, Mepilex and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2021 Mölnlycke Health Care AB. All rights reserved. HQIM00 2958

